



# 5X PCR Master Mix

Cat. ADPMX01-100  
Qty. 500 µl



## Standard Protocol

1 Set up each reaction as follows:

component	Vol./reaction	Final Conc.
5X PCR Master Mix	5 µl	1X
Primer A	Variable	0.1-1.0 µM
Primer B	Variable	0.1-1.0 µM
Distilled Water	Variable	-----
Template DNA	Variable	Variable
TOTAL Volume	25 µl	-----

2 Mix gently by pipetting the solution up and down a few times.

3 Program the thermal cycler according to the manufacturers instructions.  
For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.

4 Place the tubes in the thermal cycler and start the reaction.

## Content and Storage

Content	Cat	Product	Qty
	ADPMX01-100	5X PCR Master Mix	500 µl
Storage	stored in a stable -20 °C		

## Composition of 5X PCR Master Mix

- 150 mM Tris-HCl pH 8.5,
- 40 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>,
- 1.5 mM MgCl<sub>2</sub>,
- 0.2% Tween 20
- 0.4 mM dNTPs0
- 0.05 units/µ Taq DNA polymerase

## Description

5X PCR Master Mix is a ready-to-use 5.0x reaction mix. Simply add primers, template, and water to successfully carry out primer extensions and other molecular biology applications. Engine Taq polymerase, the NH<sub>4</sub><sup>+</sup> buffer system, dNTPs and magnesium chloride are conveniently present in the Taq DNA Polymerase Master Mix . Taq DNA Polymerase Master Mix , offers several advantages. Set up time is significantly reduced. The chance of contaminating component stocks is eliminated. Reduction of reagent handling steps leads to better reproducibility. Standard tests can be set up with the confidence that results will be consistent every time.

## NOTE

Set up reaction mixtures in an area separate from that used for DNA preparation or product analysis.

The table below shows the reaction set up for a final volume of 25 µl.

After primer extension, a sample (10 to 30% of the reaction) can be loaded directly on a gel for analysis.



## 2X PCR Dye Master Mix

Cat. ADPMX02D-100  
Qty. 1.25ml

Ver. 2024 - 07



### Standard Protocol

- 1 Set up each reaction as follows:  
Table 1 : Reaction components

Component	Vol./reaction*	Final Conc.*
2X Master Mix	25 µl	1X
Primer A (10uM)	1 µl (0.5 -5ul)	0.2 µM ( 0.1-1.0 µM)
Primer B (10uM)	1 µl (0.5 -5ul)	0.2 µM ( 0.1-1.0 µM)
Distilled Water	X µl	-----
Template DNA	X µl	Variable
25mM MgCL2	0ul (0-6ul)	1.5mM (1.5-4.5 mM)
TOTAL Volume	50 µl	-----

\* Suggest starting conditions ; theoretically used conditions brackets

- 2 Mix gently by pipetting the solution up and down a few times.
- 3 Program the thermal cycler according to the manufacturer's instructions. See table 2 for an example.  
For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.
- 4 Place the tubes in the thermal cycler and start the reaction.

### Content and Storage

Cat	Product	Qty	Store at -20 °C
ADPMX02D-100	2X PCR Dye Mix	1.25ml	

### Composition of 2X PCR Master Dye Mix

150 mM Tris-HCl pH 8.5,  
40 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>,  
1.5 mM MgCl<sub>2</sub> (final concentration)  
0.2% Tween 20  
0.4 mM dNTPs  
0.05 units/µL Taq DNA polymerase  
Inert red dye and a stabilizer

### Description

Taq DNA Polymerase 2x Master Mix RED is a ready-to-use 2x reaction mix, the NH<sub>4</sub><sup>+</sup> buffer system, dNTPs and magnesium chloride present.  
Each reaction requires 50 µl of the 2x Master Mix RED. Simply add primers, template and water to a total reaction volume of 25 µl to successfully carry out primer extensions and other molecular biology applications.  
Taq DNA Polymerase 2x Master Mix RED offers several advantages. Set up time is significantly reduced. The chance of contaminating component stocks is eliminated. Reduction of reagent handling steps leads to better reproducibility. Standard tests can be set up with the confidence that results will be consistent every time. There is no need to buy and use separate loading dyes. Simply load a portion of the reaction product onto an agarose gel for electrophoresis and subsequent visualization. The red dye front runs at 1000 – 2000 bp on a 0.5 – 1.5 % agarose gel

### Standard Protocol ( continued )

Table 2 : Three-step PCR program

Cycles	Duration of cycle	Temperature
1	2-5 minutes <sup>a</sup>	95 °C
25-35	20-30 seconds <sup>b</sup>	95 °C
	20-40 seconds <sup>c</sup>	50-65 °C
	30 seconds <sup>d</sup>	72 °C
1	5 minutes	72 °C

- Denaturation step: This step is the first regular cycling event and consists of heating the reaction to 95 °C for 20–30 seconds. It causes melting of the DNA template by disrupting the hydrogen bonds between complementary bases, yielding single-stranded DNA molecules.
- Annealing step: The reaction temperature is lowered to 50–65 °C for 20–40 seconds allowing annealing of the primers to the single stranded DNA template. Typically, the annealing temperature is about 3–5 °C below the T<sub>m</sub> (melting temperature) of the primers used.
- Extension/elongation step: Taq polymerase has its optimum activity temperature at 72 °C. At this step the DNA polymerase synthesizes a new DNA strand complementary to the DNA template strand. The extension time depends on the length of the DNA fragment to be amplified. As a rule of thumb, at its optimum temperature the DNA polymerase will polymerize a thousand bases per minute.
- Final elongation: This single step is occasionally performed at a temperature of 72 °C for 5 minutes after the last PCR cycle to ensure that any remaining single-stranded DNA is fully extended.

## NOTE

- 1 Set up reaction mixtures in an area separate from that used for DNA preparation or product analysis. The table below shows the reaction set up for a final volume of 50  $\mu$ l. After primer extension, a sample (10 to 30% of the reaction) can be loaded directly on a gel for analysis.
- 2 The final MgCl<sub>2</sub> concentration of this 2x Taq Master Mix RED is 1.5 mM. In some applications, more than 1.5 mM MgCl<sub>2</sub> is required for best results. Use 25 mM to adjust the Mg<sup>2+</sup> concentration according to table 3.

Table 3. Additional volumn (ul) of MgCl<sub>2</sub> per 50 ul reaction

Final MgCl <sub>2</sub> conc In reaction (mM)	1.5	2.0	2.5	3.0	3.5	4.0	4.5
Volume of 25 mM MgCl <sub>2</sub>	0	1	2	3	4	5	6





# 2X Hi Fi DNA Polymerase MIX

Cat. ADPMX03-100  
Qty. 1.25 ml



## Standard Protocol

### Template :

2 x Hi Fi mix is suitable for amplifying targets up to 3 kb from the following templates:

- Genomic DNA: 10–200 ng
- Plasmid DNA : 1–5 ng
- cDNA : ~100 ng starting total RNA

### Primers :

Use 0.3 μM per primer as a general starting point. For larger amounts of template (e.g., 200 ng genomic DNA), increasing the concentration up to 0.5 μM per primer may improve yield.

### Annealing Temperature :

The annealing temperature is slightly higher than with typical PCR. The optimal annealing temperature should be ~2°C lower than the T<sub>m</sub> of the primers used. A range of 58–68°C is recommended.

### Extension Time :

As little as 30 seconds per kb is suitable for most targets. Use up to 60 seconds per kb for maximum yield.

## Content and Storage

Content	Cat ADPMX03-100	Product 2X Hi Fi DNA Polymerase MIX	Qty 1.25 ml
Storage	stored in a stable -20 °C		

## Description

Hi Fi polymerase synthesizes higher yields of product from genomic DNA, cDNA, bacterial cultures. It contains a 2.5 hours half life at 96°C and is easy to amplify PCR product at G-C rich and secondary structure. 2X Hi Fi Mix is an optimized mixture containing Hi Fi polymerase, reaction buffer, dNTP and enhancer as 2-fold concentration.

2x Hi Fi mix is designed to allow the user for quick, easy preparation of reaction mixture. The 2x Hi Fi mix can be used for amplification of PCR products up to 10-15 kb and the products can be directly cloned into T-vector.

## Standard Protocol

- 1 Thaw the 2x Hi Fi mix at room temperature. Vortex the 2x Hi Fi mix and then spin it briefly in a micro centrifuge to collect the material in the bottom of the tube.
- 2 Prepare one of the following reaction mixes on ice:

Step	Volume
2x Hi Fi mix	12.5 μl
Primer1 (20 pmol)	1-2 μl
Primer2 (20 pmol)	1-2 μl
Template	1-10 μl
ddH <sub>2</sub> O	UP to 25 μl
Total	25 μl

- 3 If necessary you can scale up your volume

### Program the thermal cycler as follows:

Step	Time	Temp	Cycle
Initial denaturation	30 sec - 2 mins	95°C	1
Denature:	15 sec – 2 mins	95°C	15-30
Anneal:	15 sec – 2 mins	T <sub>m</sub> -5 °C	
Extend:	2 mins /1kbp	68-75°C	
Final extension:	1 - 10 mins 68 - 72°C		



# 2X Hot start Master Mix

Cat. ADPMX04-100  
Qty. 1.25 ml



## Standard Protocol

1 Set up each reaction as follows:

component	Vol./reaction	Final Conc.
2X Hot Start Master mix	12.5 µl	1X
Primer A	Variable	0.1-1.0 µM
Primer B	Variable	0.1-1.0 µM
Distilled Water	Variable	-----
Template DNA	Variable	Variable
TOTAL Volume	25 µl	-----

2 Mix gently by pipetting the solution up and down a few times.

3 Program the thermal cycler as follows:

Step	Temperature ( °C )	Time	Cycle
Initial denaturation	94 – 95	10 mins	1
Denaturation	94 – 95	0.2 – 2 mins	20 – 35
Annealing	50 – 68	0.2 – 2 mins	
Extension	72	1 min / kb	
Final extension	72	1 – 10 mins	1

## Content and Storage

Content	Cat	Product	Qty
	ADPMX04-100	2X Hot start Master Mix	1.25 ml
Storage	stored in a stable -20 °C		

## Description

Hot start Taq DNA Polymerase for qPCR is designed for Real-Time PCR and Hot-start PCR. A special inhibition the reaction at room temperature until after the first denaturation step. This prevents primer-dimers and other artefacts.

2X Hot Start Master mix is optimized mixture contain of Hot Start Taq enzyme, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x Hot Start master mix is designed to allow the user for quick ,easy preparation of reaction mixture. The 2x Hot Start Master mix can be amplification PCR products up to 5 kb and the products can be directly cloning into T-vector.

## Standard Protocol

After cycling, maintain the reaction at 4°C. Samples can be stored at –20°C until use.

Analyze products using standard agarose gel electrophoresis.



# 2X Hot start Dye Master Mix

Cat. ADPMX04D-100  
Qty. 1.25 ml



## Standard Protocol

1 Set up each reaction as follows:

component	Vol./reaction	Final Conc.
2X Hot Start Master mix	12.5 µl	1X
Primer A	Variable	0.1-1.0 µM
Primer B	Variable	0.1-1.0 µM
Distilled Water	Variable	-----
Template DNA	Variable	Variable
TOTAL Volume	25 µl	-----

2 Mix gently by pipetting the solution up and down a few times.

3 Program the thermal cycler as follows:

Step	Temperature ( °C )	Time	Cycle
Initial denaturation	94 – 95	10 mins	1
Denaturation	94 – 95	0.2 – 2 mins	20 – 35
Annealing	50 – 68	0.2 – 2 mins	
Extension	72	1 min / kb	
Final extension	72	1 – 10 mins	1

## Content and Storage

Content	Cat	Product	Qty
	ADPMX04D-100	2X Hot start Dye Master Mix	1.25 ml
Storage	stored in a stable -20 °C		

## Description

Hot start Taq DNA Polymerase for qPCR is designed for Real-Time PCR and Hot-start PCR. A special inhibition the reaction at room temperature until after the first denaturation step. This prevents primer-dimers and other artefacts.

2X Hot Start Dye Master mix is optimized mixture contain of Hot Start Taq enzyme, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x Hot Start dye master mix is designed to allow the user for quick ,easy preparation of reaction mixture. The 2x Hot Start Dye Master mix can be amplification PCR products up to 5 kb and the products can be directly cloning into T-vector.

## Standard Protocol

After cycling, maintain the reaction at 4°C. Samples can be stored at -20°C until use.

Analyze products using standard agarose gel electrophoresis.



## 2X Pfu DNA Polymerase MIX

Cat. ADPMX05-100  
Qty. 1.25 ml



### Standard Protocol

1 Set up each reaction as follows:

component	Vol./reaction	Final Conc.
2X Pfu DNA Polymerase mix	12.5 µl	1X
Primer A	Variable	0.1-1.0 µM
Primer B	Variable	0.1-1.0 µM
Distilled Water	Variable	-----
Template DNA	Variable	Variable
TOTAL Volume	25 µl	-----

2 Mix gently by pipetting the solution up and down a few times.

3 Program the thermal cycler according to the manufacturer's instructions.  
For maximum yield and specificity, temperatures and cycling times should be optimized for each new template target or primer pair.

4 Place the tubes in the thermal cycler and start the reaction.

### Content and Storage

Content	Cat	Product	Qty
	ADPMX05-100	2X Pfu DNA Polymerase MIX	1.25 ml

Storage stored in a stable -20 °C

### Description

2X Pfu DNA Polymerase MIX is a ready-to-use 2.0x reaction mix. Simply add primers, template, and water to successfully carry out primer extensions and other molecular biology applications. *Pfu* polymerase, the NH<sub>4</sub><sup>+</sup> buffer system, dNTPs and magnesium chloride are conveniently present in the Pfu DNA Polymerase Master Mix . offers several advantages. Set up time is significantly reduced. The chance of contaminating component stocks is eliminated. Reduction of reagent handling steps leads to better reproducibility. Standard tests can be set up with the confidence that results will be consistent every time.

### Standard Protocol

Program the thermal cycler as follows:

Step	Time	Temp	Cycle
Initial denaturation	30 sec - 2 mins	95 °C	1
Denature:	15 sec - 2 mins	95 °C	15-30
Anneal:	15 sec - 2 mins	Tm-5 °C	
Extend:	2 mins /1kbp	68-75 °C	
Final extension: 1 - 10 mins 72 °C			

After cycling, maintain the reaction at 4 °C. Samples can be stored at -20 °C until use.

### Program the thermal cycler as follows:

Set up reaction mixtures in an area separate from that used for DNA preparation or product analysis.

The table below shows the reaction set up for a final volume of 25 µl.

After primer extension, a sample (10 to 30% of the reaction) can be loaded directly on a gel for analysis.