



RNAAfter

Cat. RF0100
Qty. 100 ml



Content and Storage

RNAAfter
Cat.. RF0100



Contents

Cat. RF0100
RNAAfter100 ml

Store at RT, if precipitates form, warm at 37 °C to redissolve.

Description

RNAAfter is a non-toxic reagent for storage of various animal or plant tissues, culture cells and bacteria for RNA purification without using liquid nitrogen or -70°C freezer. The sample can be stored in RNAAfter reagent for a day at 37°C, 1 week at RT, one month at 4°C and indefinitely at -20°C. The purified RNA quality is as high as stored in liquid nitrogen.

Standard Protocol

Animal tissue: weight the tissue samples, cut the tissue samples into small pieces each with < 0.5 cm thick, add the dissected tissues into the tube with 5 volumes of RNAAfter solution. (e.g., 1 g tissue need 5 ml of RNAAfter)

Plant tissue: weight the sample, cut the tissue into small pieces and add the sample into the tube with 5 volumes of RNAAfter.

Culture cells: Spin down the cells, wash the cells with PBS buffer, resuspend the cells in PBS buffer, add 5 volumes of RNAAfter.

White blood cell: Separate white blood cells from whole blood, wash with PBS buffer, resuspend the cells in PBS buffer, add 5 volumes of RNAAfter. Do not use whole blood for storage, which will precipitate the RNAAfter.

Bacteria: Spin down the cells, wash the cells with TE buffer, resuspend the cells in TE buffer, add 5 volumes of RNAAfter.

Sample Storage:
After submersing in RNAAfter, the sample can be stored for a day at 37°C, 1 week at RT, one month at 4°C and indefinitely at -20°C. Sample can be thawed and frozen many times without affecting the RNA quality. It may form crystal in lower temperature, but it will not affect the RNA purification.

Standard Protocol (continued)

RNA Purification from samples in RNAAfter:

Tissue: remove the RNAAfter solution or use forceps to take out the tissue from solution, add RNA extraction lysis solution(e.g., Azol RNA purification reagent,Cat : Azol.200) , proceed to standard protocols.

Cells: Spin down the cell at 5000 x g for 3 min, remove the RNAAfter, add RNA extraction lysis solution (e.g., Azol RNA purification reagent,Cat : Azol.200), proceed to standards protocols.