


ARROW Prestained Protein Ladder



Content and Storage

Content	Cat	Product	Qty
	APM0671.250	Prestained Protein Ladder	250 µl /vial
Storage	Stable for up to 2 weeks at 25°C Stable for up to 3 months at 4°C For long term storage, store at -20°C		

Description

The **ARROW Prestain Protein Marker** is designed for visualization of protein molecular of Western transfer efficiency. The Prestained Ladder consists of 10 prestained proteins Covering the range 10-180kDa.

Proteins are covalently coupled to A blue dye, except the fourth protein band from the top is coupled To a green dye, the eighth protein band from the top is coupled A red dye.

Quality Control

5 µl of ARROW Prestained Protein Ladder resolves 10 bands in 4-20% SDS-PAGE (Tris-glycine buffer) and after Western blotting to PVDF membrane.